


VHP Flexible Film Isolator System

(Positive or Negative Pressure)

Isolators can be operating in hours vs. days!



Pictured above are two germ-free, flexible film, positive pressure isolator systems equipped for vapor phase hydrogen peroxide (VHP) sterilization.

Specially Equipped for Vapor Phase Hydrogen Peroxide Sterilization

(With VHP sterilization, isolators can be up and operating in hours instead of days.)

VHP Isolator System features:

- 1 1/2" steel powder coated single shelf frame to prevent bending or warping. (optional adjustable frame with motors.)
- HDPE tabletop with inserts for attaching isolator components.
- Two-way locking wheels.
- 5/8" PVC hold-down rod and SS clamps.
- Blower and manifold system.
- (1) Durable 20 mil vinyl isolator.
- (1) 20 mil urethane floor mat protectors for inside isolators.
- 3/8" closed cell foam pads for under isolators.
- (1) complete 5" diameter inlet filter & (1) complete 5" diameter inlet filter. Includes: filter media and filter covers, rubber bands and clamps.
- (1) pair of 12" diameter/8" CBC glove port extensions for better worker dexterity and less stress on isolator. (Note: Available without.)
- (1) pairs of 15 mil Hypalon/CSM dry box gloves with glove ring kit for attaching gloves to isolator (other types available).
- (1) Transfer port. (size and type to your specifications.)
- (1) Port to isolator neoprene band and clamp for sealing isolator to transfer port.
- (2) 1" neoprene stoppers (other types available).
- (2) internal racks that can accommodate up to 9 mouse cages each.
- (1) outside transfer port cap.
- (1) inside transfer port cap.
- (1) patent clamp.

For more information about our CBC germ-free, flexible film, VHP isolator systems or other equipment used for gnotobiotic, transgenic or microbiome research, call +49 7681 4979 395 or e-mail: av@cbclean.com



Quality Built to Protect
Both Research and Researchers

Szymaniak Vits UG • Dettenbach 9 • 79183 Waldkirch, Germany